


Product Datasheet

Human CCL21/6Ckine ELISA Kit (DIY Antibody Pairs) (orb1098164)

Description	Human CCL21/6Ckine ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S24-P134
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Uniprot ID	000585
Expiration Date	Please enquire.



Human CCL21/6Ckine EZ-Set ELISA Kit stan...