


Product Datasheet

Human IGFBP3 ELISA Kit (DIY Antibody Pairs) (orb1098166)

| | |
|----------------------------|---|
| Description | Human IGFBP3 ELISA Kit (DIY Antibody Pairs) |
| Reactivity | Human |
| Range | 156 pg/ml - 10,000 pg/ml |
| Tested Applications | ELISA |
| Immunogen | Expression system for standard: NS0; Immunogen sequence: A28Q-291 |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Sample Types | cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA). |
| Uniprot ID | P17936 |
| Expiration Date | Please enquire. |



Human IGFBP3 EZ-Set ELISA Kit standard C...