


Product Datasheet

Mouse IL-1RA ELISA Kit (DIY Antibody Pairs) (orb1098168)

Description	Mouse IL-1RA ELISA Kit (DIY Antibody Pairs)
Reactivity	Mouse
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: R27-Q178
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Uniprot ID	P25085
Expiration Date	Please enquire.



Mouse IL-1RA EZ-Set ELISA Kit standard C...