


Product Datasheet

Rat Galectin-1 ELISA Kit (DIY Antibody Pairs) (orb1098173)

Description	Rat Galectin-1 ELISA Kit (DIY Antibody Pairs)
Reactivity	Rat
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A2-E135
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	P11762
Expiration Date	Please enquire.



Rat
Galectin-1
EZ-Set
ELISA Kit
standard...