


Product Datasheet

Human Angiopoietin-2 ELISA Kit (DIY Antibody Pairs) (orb1098179)

Description	Human Angiopoietin-2 ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Y19-F496
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva and human milk.
Uniprot ID	015123
Expiration Date	Please enquire.



Human Angiopoietin-2 EZ-Set ELISA Kit st...