


Product Datasheet

Human EPCR ELISA Kit (DIY Antibody Pairs) (orb1098180)

Description	Human EPCR ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S18-S210
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Uniprot ID	Q9UNN8
Expiration Date	Please enquire.



HUMAN
EPCR EZ-
Set ELISA
Kit
standard
cur...