

Product Datasheet

Human Angiopoietin-1 ELISA Kit (DIY Antibody Pairs) (orb1098182)

Description Human Angiopoietin-1 ELISA Kit (DIY Antibody Pairs)

Reactivity Human

Range 156 pg/ml - 10,000 pg/ml

Tested Applications ELISA

Immunogen Expression system for standard: NS0; Immunogen sequence: H16-F498


Storage Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)

Note For research use only

Sample Types cell culture supernatants, serum, plasma (heparin, EDTA), saliva and human milk.

Uniprot ID [Q15389](#)

Expiration Date Please enquire.



Human Angiopoietin-1 EZ Set ELISA Kit st...