


Product Datasheet

Human MIF ELISA Kit (DIY Antibody Pairs) (orb1098185)

Description	Human MIF ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: P2-A115
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and human milk.
Uniprot ID	P14174
Expiration Date	Please enquire.



Human MIF
EZ-Set
ELISA Kit
standard
curv...