


## Product Datasheet

# Human IL-1 beta ELISA Kit (DIY Antibody Pairs) (orb1098186)

<b>Description</b>	Human IL-1 beta ELISA Kit (DIY Antibody Pairs)
<b>Reactivity</b>	Human
<b>Range</b>	3.9 pg/ml - 250 pg/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: E.coli; Immunogen sequence: A117-S269
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Sample Types</b>	cell culture supernatants, serum and plasma (heparin, EDTA, citrate)
<b>Uniprot ID</b>	<b>P01584</b>
<b>Expiration Date</b>	Please enquire.



Human IL-1  
beta EZ-  
Set ELISA  
Kit  
standar...