


Product Datasheet

Human KLK1 ELISA Kit (DIY Antibody Pairs) (orb1098193)

Description	Human KLK1 ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: I25-S262
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Uniprot ID	P06870
Expiration Date	Please enquire.



Human
KLK1 EZ-
Set ELISA
Kit
standard
cur...