


Product Datasheet

Rat HE4/WFDC2 ELISA Kit (DIY Antibody Pairs) (orb1098199)

Description	Rat HE4/WFDC2 ELISA Kit (DIY Antibody Pairs)
Reactivity	Rat
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: E31-F168
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q8CHN3
Expiration Date	Please enquire.



Rat
HE4/WFDC2
EZ-Set
ELISA Kit
standard ...