


Product Datasheet

Human Diablo/SMAC ELISA Kit (DIY Antibody Pairs) (orb1098203)

| | |
|----------------------------|---|
| Description | Human Diablo/SMAC ELISA Kit (DIY Antibody Pairs) |
| Reactivity | Human |
| Range | 31.2 pg/ml - 2,000 pg/ml |
| Tested Applications | ELISA |
| Immunogen | Expression system for standard: E.coli; Immunogen sequence: A56-D239 |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Sample Types | cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA). |
| Uniprot ID | Q9NR28 |
| Expiration Date | Please enquire. |



Human
Diablo/SMAC
EZ-Set
ELISA Kit
stand...