


Product Datasheet

Rat Laminin ELISA Kit (DIY Antibody Pairs) (orb1098208)

Description	Rat Laminin ELISA Kit (DIY Antibody Pairs)
Reactivity	Rat
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: from rat sarcoma basement membrane
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Uniprot ID	D4A409, F1MAA7, A0A0G2K023, F1LRH4, F1M614, F1MAN8, D3ZQN7, F1LPI5, F1LTF8, D3ZN05, P15800
Expiration Date	Please enquire.



Rat Laminin EZ-Set ELISA Kit standard cu...