


Product Datasheet

Human IL-1RA ELISA Kit (DIY Antibody Pairs) (orb1098215)

Description	Human IL-1RA ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: R26-E177
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Uniprot ID	P18510
Expiration Date	Please enquire.



Human IL-1RA EZ-Set ELISA Kit standard C...