


Product Datasheet

Mouse IFN gamma ELISA Kit (DIY Antibody Pairs) (orb1098226)

Description	Mouse IFN gamma ELISA Kit (DIY Antibody Pairs)
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: H23-C155
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum
Uniprot ID	P01580
Expiration Date	Please enquire.



Mouse IFN
gamma
EZ-Set
ELISA Kit
standar...