


Product Datasheet

Rat MIP-1Alpha/CCL3 ELISA Kit (DIY Antibody Pairs) (orb1098230)

Description	Rat MIP-1Alpha/CCL3 ELISA Kit (DIY Antibody Pairs)
Reactivity	Rat
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A24-A92
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	10335 MW
Uniprot ID	P50229
Expiration Date	Please enquire.



Rat MIP-1alpha/CCL3 EZ Set ELISA Kit sta...