


Product Datasheet

Human Tissue Factor/F3 ELISA Kit (DIY Antibody Pairs) (orb1098234)

Description	Human Tissue Factor/F3 ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S33-E251
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA, citrate) , Cell lysates, and urine.
MW	33068 MW
Uniprot ID	P13726
Expiration Date	Please enquire.



Human
Tissue
factor/F3
EZ Set
ELISA Kit
...