


Product Datasheet

Human ADAMTS13 ELISA Kit (DIY Antibody Pairs) (orb1098235)

Description	Human ADAMTS13 ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: Q34-W688
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	153604 MW
Uniprot ID	Q76LX8
Expiration Date	Please enquire.



Human
ADAMTS13
EZ Set
ELISA Kit
standard...