


Product Datasheet

Human CD40/TNFRSF5 ELISA Kit (DIY Antibody Pairs) (orb1098244)

| | |
|----------------------------|---|
| Description | Human CD40/TNFRSF5 ELISA Kit (DIY Antibody Pairs) |
| Reactivity | Human |
| Range | 15.6 pg/ml - 1,000 pg/ml |
| Tested Applications | ELISA |
| Immunogen | Expression system for standard: NS0; Immunogen sequence: E21-R193 |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Sample Types | cell culture supernatants, serum, plasma (heparin, EDTA) , saliva, urine and human milk. |
| MW | 30619 MW |
| Uniprot ID | P25942 |
| Expiration Date | Please enquire. |



Human CD40/TNFRSF5 EZ Set ELISA Kit stan...