


## Product Datasheet

# Human Angiogenin/ANG ELISA Kit (DIY Antibody Pairs) (orb1098251)

<b>Description</b>	Human Angiogenin/ANG ELISA Kit (DIY Antibody Pairs)
<b>Reactivity</b>	Human
<b>Range</b>	78 pg/ml - 5,000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: E.coli; Immunogen sequence: Q25-P147
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Sample Types</b>	cell culture supernatants, cell lysates, tissue homogenates, serum and plasma (heparin, EDTA).
<b>MW</b>	16550 MW
<b>Uniprot ID</b>	<a href="#">P03950</a>
<b>Expiration Date</b>	Please enquire.



Human Angiogenin/ANG EZ Set ELISA Kit st...