

Product Datasheet

Human Lipocalin-2/NGAL ELISA Kit (DIY Antibody Pairs) (orb1098252)


Description	Human Lipocalin-2/NGAL ELISA Kit (DIY Antibody Pairs) (antibody pairs and standards for assay development). Quantitate Human LCN2 in cell culture supernatants, serum, plasma (heparin), saliva and urine.
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q21-G198
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin), saliva and urine.
MW	22588 Da
Uniprot ID	P80188
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



Human Lipocalin-2/NGAL EZ Set ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)