

Product Datasheet

SARS-CoV-2 S1 protein Antibody (orb1100528)

Description	anti- 2019 nCoV S1 protein antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	S1-RBD
Target	SARS-CoV-2 S1 protein
Form/Appearance	liquid
Storage	PBS pH 7.4,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	P0DTC2
Dilution Range	ELISA: Use at an assay dependent concentration(Suggestions: Use a concentration of 0.1 - 1µg/ml);WB: 1:1000-1:10000
Expiration Date	12 months from date of receipt.