

Product Datasheet

anti- NFATC3 antibody (orb1100597)

Description	anti- NFATC3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 3
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	130-170 kDa
Uniprot ID	Q12968
Dilution Range	WB: 1:200-1:2000; IHC: 1:50-1:200
Expiration Date	12 months from date of receipt.