

Product Datasheet

anti- TAB3 antibody (orb1100616)

Description	anti- TAB3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	TGF-beta-activated kinase 1 and MAP3K7-binding protein 3
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	75-100 kDa
Uniprot ID	Q8N5C8
Dilution Range	WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:20-1:100
Expiration Date	12 months from date of receipt.