


Product Datasheet

Rat Gpnmb (Transmembrane glycoprotein NMB) ELISA Kit (orb1100751)

Description	Rat Gpnmb (Transmembrane glycoprotein NMB) ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/ml
Target	Gpnmb
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q6P7C7
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Rat Gpnmb
(Transmembrane
glycoprotein
NM...