


## Product Datasheet

# Rat Cd200 (OX-2 membrane glycoprotein) ELISA Kit (orb1100757)

<b>Description</b>	Rat Cd200 (OX-2 membrane glycoprotein) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	78.125-5000pg/ml
<b>Target</b>	Cd200
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>P04218</b>
<b>Sensitivity</b>	46.875pg/ml
<b>Expiration Date</b>	Please enquire.



Rat Cd200 (OX-2 membrane glycoprotein) E...