


Product Datasheet

Rat ALB (urine albumin) ELISA Kit (orb1100771)

Description	Rat ALB (urine albumin) ELISA Kit
Reactivity	Rat
Range	7.813-500ng/ml
Target	ALB
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Sensitivity	4.688ng/ml
Expiration Date	Please enquire.



Rat ALB
(urine
albumin)
ELISA Kit