


Product Datasheet

Rat CD200 R1 (OX-2 membrane glycoprotein receptor) ELISA Kit (orb1100786)

| | |
|------------------------|--|
| Description | Rat CD200 R1 (OX-2 membrane glycoprotein receptor) ELISA Kit |
| Reactivity | Rat |
| Range | 31.25-2000pg/ml |
| Target | CD200 R1 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Sample Types | Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples |
| Uniprot ID | Q9ES58 |
| Sensitivity | 18.75pg/ml |
| Expiration Date | Please enquire. |



Rat CD200 R1 (OX-2 membrane glycoprotein...)