


Product Datasheet

Rat UBL4A (Ubiquitin-like protein 4A) ELISA Kit (orb1100791)

Description	Rat UBL4A (Ubiquitin-like protein 4A) ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	UBL4A
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	B2GV38
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat UBL4A (Ubiquitin-like protein 4A) EL...