


Product Datasheet

Rat MMP14 (Matrix metalloproteinase-14) ELISA Kit (orb1100799)

Description	Rat MMP14 (Matrix metalloproteinase-14) ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	MMP-14
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q10739
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat MMP14
(Matrix
metalloproteinase-
14) ...