


## Product Datasheet

# Rat MLCK/MYLK (Myosin light chain kinase) ELISA Kit (orb1100802)

<b>Description</b>	Rat MLCK/MYLK (Myosin light chain kinase) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	31.25-2000pg/ml
<b>Target</b>	MLCK/MYLK
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>P20689</b>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



Rat  
MLCK/MYLK  
(Myosin  
light chain  
kinase...