


Product Datasheet

Mouse Aplnr (Apelin receptor) ELISA Kit (orb1100822)

Description	Mouse Aplnr (Apelin receptor) ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/ml
Target	Aplnr
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q9WV08
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Mouse
Aplnr
(Apelin
receptor)
ELISA Kit