


## Product Datasheet

# Mouse UALB (Urinary Albumin)ELISA Kit (orb1100827)

<b>Description</b>	Mouse UALB (Urinary Albumin)ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.781-50ng/ml
<b>Target</b>	UALB
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Sensitivity</b>	0.469ng/ml
<b>Expiration Date</b>	Please enquire.



Mouse  
UALB  
(Urinary  
Albumin)  
ELISA Kit