


Product Datasheet

Mouse AT1R (angiotensin II receptor type 1) ELISA Kit (orb1100868)

Description	Mouse AT1R (angiotensin II receptor type 1) ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/ml
Target	AT1R
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q9WVK0
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Mouse AT1R
(angiotensin
II receptor
typ...