


Product Datasheet

Mouse ACE2 (Angiotensin I Converting Enzyme 2) ELISA Kit (orb1100870)

Description	Mouse ACE2 (Angiotensin I Converting Enzyme 2) ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/ml
Target	ACE2
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q8R0I0
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Mouse ACE2
(Angiotensin
I Converting
Enz...