


Product Datasheet

Human IGSF11(Immunoglobulin superfamily member 11) ELISA Kit (orb1100896)

Description	Human IGSF11(Immunoglobulin superfamily member 11) ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	IGSF11
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q5DX21
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Human IGSF11(Immunoglobulin superfamily ...