


Product Datasheet

Human FIS1 (Mitochondrial fission 1 protein) ELISA Kit (orb1100913)

| | |
|------------------------|--|
| Description | Human FIS1 (Mitochondrial fission 1 protein) ELISA Kit |
| Reactivity | Human |
| Range | 0.156-10ng/ml |
| Target | FIS1 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | Q9Y3D6 |
| Sensitivity | 0.094ng/ml |
| Expiration Date | Please enquire. |



Human FIS1
(Mitochondrial
fission 1
prot...