


## Product Datasheet

# Human 10-HETE(10-Hydroxyeicosatetraenoic Acid) ELISA Kit (orb1100972)

<b>Description</b>	Human 10-HETE(10-Hydroxyeicosatetraenoic Acid) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/ml
<b>Target</b>	10-HETE
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sensitivity</b>	0.094ng/ml
<b>Expiration Date</b>	Please enquire.



Human 10-HETE(10-Hydroxyeicosatetraenoic...