


Product Datasheet

Human RHEBL1 (GTPase RhebL1) ELISA Kit (orb1100997)

Description	Human RHEBL1 (GTPase RhebL1) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	RHEBL1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q8TAI7
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
RHEBL1
(GTPase
RhebL1)
ELISA Kit