


Product Datasheet

Human DNM1L (Dynamin-1-like protein) ELISA Kit (orb1101011)

Description	Human DNM1L (Dynamin-1-like protein) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	DNM1L
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O00429
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
DNM1L
(Dynamin-
1-like
protein)
ELI...