


Product Datasheet

Human HTR2A (5-hydroxytryptamine receptor 2A) ELISA Kit (orb1101018)

| | |
|------------------------|---|
| Description | Human HTR2A (5-hydroxytryptamine receptor 2A) ELISA Kit |
| Reactivity | Human |
| Range | 31.25-2000pg/ml |
| Target | HTR2A |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Uniprot ID | P28223 |
| Sensitivity | 18.75pg/ml |
| Expiration Date | Please enquire. |



Human HTR2A (5-hydroxytryptamine recepto...