


Product Datasheet

Human TSPAN12(Tetraspanin-12) ELISA Kit (orb1101026)

Description	Human TSPAN12(Tetraspanin-12) ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	TSPAN12
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	095859
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Human
TSPAN12(Tetraspanin-
12) ELISA Kit