


Product Datasheet

Human HS-IL-17(High sensitive Interleukin 6) Accquant ELISA Kit (orb1101087)

Description	Human HS-IL-17(High sensitive Interleukin 6) Accquant ELISA Kit
Reactivity	Human
Range	6.25-400pg/ml
Target	IL-17
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q16552
Sensitivity	3.75pg/ml
Expiration Date	Please enquire.



Human HS-IL-17(High sensitive Interleuki...