


Product Datasheet

Human HS-IL-18(High sensitive Interleukin 18) Accquant ELISA Kit (orb1101095)

Description	Human HS-IL-18(High sensitive Interleukin 18) Accquant ELISA Kit
Reactivity	Human
Range	1.563-100pg/ml
Target	IL-18
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q14116
Sensitivity	0.983pg/ml
Expiration Date	Please enquire.



Human HS-IL-18(High sensitive Interleuki...