


## Product Datasheet

### m04 (MCMV) antibody (orb1101105)

<b>Description</b>	Mouse monoclonal antibody to m04 (MCMV)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FACS, IF, IHC, IP, PLA
<b>Immunogen</b>	recombinant m04-Fc fusion protein expressed in HEK293T cells (extracellular domain without signal sequence - amino-acids at the positions 24-219)
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 mg/mL
<b>Storage</b>	Long term -20 °C, short term +4 °C. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Affinity purified
<b>Expiration Date</b>	12 months from date of receipt.



Proximity ligation assay demonstrating i...