


Product Datasheet

M112-113/E1 (MCMV) antibody (orb1101109)

Description	Mouse monoclonal antibody to M112-113/E1 (MCMV)
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	The immunogen consisted of MCMV virions
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Long term -20 °C, short term +4 °C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Affinity purified
Expiration Date	12 months from date of receipt.



WB analysis
of cells:
a=uninfected
cells...