


Product Datasheet

m152 (MCMV) antibody (orb1101112)

Description	Mouse monoclonal antibody to m152 (MCMV)
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, WB
Immunogen	The immunogen consisted of N-terminal portion of m152
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Long term -20 °C, short term +4 °C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG2c
Clonality	Monoclonal
Purity	Affinity purified
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis
ofMEF-K
cells u...