


Product Datasheet

M55/gB (MCMV) antibody (orb1101116)

Description	Mouse monoclonal antibody to M55/gB (MCMV)
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	FACS, WB
Immunogen	MCMV virions
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Long term -20 °C, short term +4 °C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Affinity purified
Expiration Date	12 months from date of receipt.



WB analysis
of cells:
a=uninfected
cells...