


Product Datasheet

MATp1 (MCMV) antibody (orb1101119)

Description	Mouse monoclonal antibody to MATp1 (MCMV)
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IP, PLA, WB
Immunogen	His-tagged MATp1 produced in E.coli
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Long term -20 °C, short term +4 °C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Affinity purified
Expiration Date	12 months from date of receipt.



WB analysis of uninfected cells, WT MCMV...