


Product Datasheet

ORF31/gB (VZV) antibody (orb1101139)

Description	Mouse monoclonal antibody to ORF31/gB (VZV)
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The immunogen consisted of ORF31 C terminal domain (amino acids 750-868) and was produced in E.coli
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Long term -20 °C, short term +4 °C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Affinity purified
Expiration Date	12 months from date of receipt.



WB analysis of cells:
a=uninfected cells...